Centrius



Low Profile - Alternating Mattress Replacement

The Intelligent Pressure Sensing Technology (IPS), responds to patient movements on the mattress by automatically adjusting internal mattress pressure, this allows Centrius to regulate airflow and continuously provide total envelopment.





Alternating Pressure Relief



CLP Therapy Pressure Redistribution



1300 547 877 | keystonehealth.care

Keystone Healthcare

MAX Inflation



Whisper Quiet



Lightweight Power Unit



Memory Recall



1 Power/Mute/Lock	
2 Panel Lock	System is locked when this light is orange

3 Comfort Level Settings Auto adjust to individual patient comfort

4 Mode Button Changes the therapy mode

5 CLP Therapy All air cells maintain constant pressure support

6 Max Inflate Rapidly hyper-inflates the mattress providing a firm and convenient surface for nursing tasks

7 Dynamic Mode 2:1 alternating therapy, 10 minute cycle

8 Low Pressure Alert Orange LED indicates insufficient pressure

Intelligent Pressure Sensing Responds to patient movements on mattress by automatically

adjusting internal mattress pressure

Power Failure Alarm

Audible and visual power failure alarms to alert caregivers

in case of abnormal power outage

Hands Free Automatic System

Soft, Medium or Firm fully automated comfort modes, with

auto upright support

Transport Mode Connect the two quick couplings together to retain air in

the mattress, can stay inflated for up to 15 hours

TECHNICAL DETAILS

GENERAL FEATURES

Risk/Category High to Very High
Waterlow Score 15+ to 20+

Pressure Sore Grading 2 to 4

Safe Working Load 16 to 250kg

Compatible Support Surfaces 5" Overlay, 7", Lux and FIT Mattress Replacements -

Single and King Single

Travel/Weight Pump weighs 2.3kg

Clinical Settings Suitable for Hospital, Home Care and Nursing Home environments

Warranty 2 year

All Carilex Medical ANZ products are registered with the Therapeutic Goods Administration.

