## Envelop



## Modular Mattress Replacement





**Mattress Replacement** 



Cable Management



Infection Control



Rapid CPR Deflation



Quick Connect Coupling



**Durable Material** 



8 inch Square Air Cells



High MVP Rate



**Easy Cleaning** 



## **ENVELOP**

## Modular Mattress Replacement

The Carilex Mattress Replacements are designed to facilitate blood circulation and provide maximum tissue pressure redistribution and relief for treating and preventing pressure injuries.

They form the ideal therapy system when used with a compatible pump, providing low noise levels and a comfortable modular cell design for superior sleeping quality.

GENERAL FEATURES	
S31 tested	79% immersion, creating the best envelopment effect in the industry
Rapid CPR Pull Tab	Deflating the entire mattress within seconds in case of an emergency
Bariatric Extension	Side extensions fitting width of 90cm, 105cm and 120cm – one mattress fits all
Square Air Cells	No gap between air cells delivers maximum support and pressure redistribution while patient is lying down. This also gives additional safety for egress and ingress.
Directional Snap	Specifically designed buttons provide better structural stability from air cells
Multiple Coverlet Material	Standard breathable quilted nylon coverlet or optional Dartex coverlet providing high quality waterproof and air permeable medical fabric
Transportation Plug	Allows the mattress to stay inflated after being detached from the power unit
TECHNICAL DETAILS	
Risk/Category	V 18.1
mak/ category	Very High
Waterlow Score	very High 20+
5 ,	
Waterlow Score	20+
Waterlow Score Pressure Sore Grading	20+ 4
Waterlow Score Pressure Sore Grading Safe Working Load	20+ 4 Up to 450kg
Waterlow Score Pressure Sore Grading Safe Working Load Air Cell Configuration	20+ 4 Up to 450kg 20 x 8" square cells
Waterlow Score Pressure Sore Grading Safe Working Load Air Cell Configuration Cycle Times	20+  4  Up to 450kg  20 x 8" square cells  10 minutes, 2:1  Single 90cm x 200cm x 25cm
Waterlow Score Pressure Sore Grading Safe Working Load Air Cell Configuration Cycle Times Available Size/Dimensions	20+  4  Up to 450kg  20 x 8" square cells  10 minutes, 2:1  Single 90cm x 200cm x 25cm King Single 107cm x 200cm x 25cm
Waterlow Score Pressure Sore Grading Safe Working Load Air Cell Configuration Cycle Times Available Size/Dimensions Securing Mattress to Bed	20+  4  Up to 450kg  20 x 8" square cells  10 minutes, 2:1  Single 90cm x 200cm x 25cm King Single 107cm x 200cm x 25cm Buckle Straps
Waterlow Score Pressure Sore Grading Safe Working Load Air Cell Configuration Cycle Times Available Size/Dimensions Securing Mattress to Bed Travel Weight	20+  4  Up to 450kg  20 x 8" square cells  10 minutes, 2:1  Single 90cm x 200cm x 25cm King Single 107cm x 200cm x 25cm  Buckle Straps  Mattress rolls up into a bag, weighing 11.5kg

Manufacturing for over 30 Years - 10 Years in the Australian Market

