

Dual Plus

Modular Mattress Replacement





Mattress Replacement



Alternating Pressure Relief



8 inch Air Cells



Patient Envelopment



Pressure Redistribution



DTX Coverlet



TPU Soft Air Cells



Quick Connect Coupling



Easy Cleaning



DUAL PLUS Modular Mattress Replacement

The Dual Plus mattress is customisable to suit individual patient needs. Localised healing can be created by easily swapping any of the 8" cells with the 5" bridging cell, providing reduced pressure zones, for example, heel zone.

GENERAL FEATURES	
Easily Exchangeable Air Cells	All cells are connected with CPC connectors and are easy to repair on site
Instant CPR Dial	Rapid deflating the entire mattress within seconds in case of an emergency
Extra Soft 100% TPU Material	Ensuring patients safety from sharp and protruding objects while providing optimal comfort
Excellent Hygiene Control	Durable and washable coverlet that completely covers the mattress. Hospital grade mattress base with easily wipe clean material
Directional Snap	Specifically designed buttons, provide better structural stability from air cells
Easy Carry Design	Mattress is designed specifically for easy carrying
Bridging Cells	Green 5" bridging cell reduce interface pressure and are used to create localised healing zones. 1 per mattress, additional bridging cells can be purchased
TECHNICAL DETAILS	
Risk/Category	Very High
	Very High 20+
Risk/Category	, •
Risk/Category Waterlow Score	20+
Risk/Category Waterlow Score Pressure Sore Grading	20+ 4
Risk/Category Waterlow Score Pressure Sore Grading Safe Working Load	20+ 4 Up to 250kg
Risk/Category Waterlow Score Pressure Sore Grading Safe Working Load Air Cell Configuration	20+ 4 Up to 250kg 20 x 8" cells
Risk/Category Waterlow Score Pressure Sore Grading Safe Working Load Air Cell Configuration Cycle Times	20+ 4 Up to 250kg 20 x 8" cells 10 minutes, 2:1
Risk/Category Waterlow Score Pressure Sore Grading Safe Working Load Air Cell Configuration Cycle Times Available Size/Dimensions	20+ 4 Up to 250kg 20 x 8" cells 10 minutes, 2:1 Single 90cm x 200cm x 21cm, King Single 107cm x 200cm x 21cm
Risk/Category Waterlow Score Pressure Sore Grading Safe Working Load Air Cell Configuration Cycle Times Available Size/Dimensions Securing Mattress to Bed	20+ 4 Up to 250kg 20 x 8" cells 10 minutes, 2:1 Single 90cm x 200cm x 21cm, King Single 107cm x 200cm x 21cm Buckle Straps
Risk/Category Waterlow Score Pressure Sore Grading Safe Working Load Air Cell Configuration Cycle Times Available Size/Dimensions Securing Mattress to Bed Travel/Weight	20+ 4 Up to 250kg 20 x 8" cells 10 minutes, 2:1 Single 90cm x 200cm x 21cm, King Single 107cm x 200cm x 21cm Buckle Straps Mattress rolls up into a bag, weighing 11.3kg

All Carilex Medical ANZ products are registered with the Therapeutic Goods Administration.

